

Overbed Tables

Model 149-BAR




Features:

- Bariatric Overbed Table
- High pressure laminate 18" x 40" top
- Structurally reinforced with a tubular steel welded frame which extends virtually the entire length of the table top
- Base: Spring assisted height mechanism with strong positive lock in 1" increments between 30.5" and 45.5"
- H style base design
- 35" opening to use with most bariatric chairs Structurally reinforced to meet the higher strength requirements
- Standard Opal powder coated column & base
- Standard Gunstock Walnut Finish

Specifications:

- Top: 18" x 40"
- Height: 30.5" - 45.5"
- 2" Casters
- Weight capacity: 500 lbs evenly distributed
- Warranty: 1 Year



Finish Options



Bannister
Oak

Clear Maple

English Oak

Figured
Mahogany

Fusion
Maple



Gunstock
Walnut
(Standard)


Nepal Teak

Quartered
Oak

Solar Oak

Wild Cherry

Standard Base Finish Choices



Brown

Gray

Opal
(Standard)

Black

Chrome

*Additional costs may apply to finish upgrades