

## Overbed Tables

### *Model 149-BAR*




### Features:

- Bariatric Overbed Table
- High pressure laminate 18" x 40" top
- Structurally reinforced with a tubular steel welded frame which extends virtually the entire length of the table top
- Base: Spring assisted height mechanism with strong positive lock in 1" increments between 30.5" and 45.5"
- H style base design
- 35" opening to use with most bariatric chairs
- Structurally reinforced to meet the higher strength requirements
- Standard Opal powder coated column & base
- Standard Gunstock Walnut Finish

### Specifications:

- Top: 18" x 40"
- Height: 30.5" - 45.5"
- 2" Casters
- Weight capacity: 500 lbs evenly distributed
- Warranty: 1 Year



## Finish Options



Bannister  
Oak

Clear Maple

English Oak

Figured  
Mahogany

Fusion  
Maple



Gunstock  
Walnut


Nepal Teak

Quartered  
Oak

Solar Oak

Wild Cherry

### Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome

\*Additional costs may apply to finish upgrades