



Adult Beds

Model CC110

Novum Medical Products, Inc. combines convenience, aesthetics and quality in a well-built line of adult beds. Our beds have many features found in costlier beds at an economic price. Contemporary style head and foot boards will complement any decor.

Specifications:

Mattress Deck Length: 76" / 80"

*optional 4" extender available

Mattress Deck Width: 35"

*optional 42" extender available

Overall Bed Length: 85.5"

Overall Bed Width: 39.75"

Maximum Back Angle: 70°

Maximum Knee Angle: 30°

Minimum Height: 8.9"

Horizontal Translation: less than 1"

Safe Workload: 500 lbs.

