



**NOVUM**  
Medical Products

# C100 Mattress



## Features:

- Comes standard in blue Rec 5 vinyl top with birthday bear vinyl bottom
- 2" of 1820 softer polyfoam on top of 2" of 1845 standard foam
- Waterproof cover
- Perforated to work well with gatch deck to allow head and knee bend
- Latex free
- Does not contain any plastics or plasticizers
- Antistatic
- Non-allergenic fabrics
- CAL 117 & CAL 129 compliant
- Warranty: 1 year

## Options:

- Also available in Baby-soft white
- Can be made in other sizes upon request

Phone | 716.759.7200  
Toll Free | 800.274.2742  
Fax | 716.759.0700  
Email | [info@novummed.com](mailto:info@novummed.com)  
[www.novummed.com](http://www.novummed.com)