



C500 Mattresses



Features:

- Standard 1845 foam
- Comes standard in Birthday Bear vinyl print
- Also available in Baby-soft white
- Waterproof cover
- Latex free
- Does not contain any plastics or plasticizers
- Antistatic
- Non-allergenic fabrics
- CAL 117 & CAL 129 compliant
- Warranty: 1 year

Model #	Dimensions	Description
C503	36" (L) x 24" (W) x 3" (H)	Foam Mattress For Neonatal Crib
C505	44" (L) x 30" (W) x 3" (H)	Foam Mattress For Infant Crib
C506	44" (L) x 30" (W) x 4" (H)	Foam Mattress For Infant Crib
C507	60" (L) x 30" (W) x 3" (H)	Foam Mattress For Child Crib
C508	60" (L) x 30" (W) x 4" (H)	Foam Mattress For Child Crib
C509	72" (L) x 36" (W) x 3" (H)	Foam Mattress For Youth Crib
C510	72" (L) x 36" (W) x 4" (H)	Foam Mattress For Youth Crib