



Cage Crib Top

Best used for very active patients who may be at risk to climb out of a standard crib. The cage crib top helps to keep the child safely in the crib.

Features:

- Powder-coated steel, confetti finish
- Sides slide up by squeezing two black tabs and pushing upwards to lock in place
- Can be mounted to any of Novum Medical's Standard and ICU cribs
- Halfway and full locking points when lifting up

Models/Specifications:

- C2000 - For 30" x 60" Child Crib.
- C2002 - For 36" x 72" Youth Crib

