

Cribettes

The 33" deck length is comfortable for individual neonates or twins and has adequate space for monitoring and observation. The Neonatal Crib is larger than a bassinet but smaller than a standard crib for easy access. The height and half rails make it more comfortable for the caregiver resulting in reduced back strain.

Standard Features:

- Easy Access Safety Gates
- 5" Casters (1 steering, 1 swivel, 2 braking)
- Powder coated confetti finish
- Solid Pan Sleeping System (non-helical design) with Fowler
- 3 Increments of Rail Stopping Points
- Crib Height: 48"
- Deck Height (from floor): 35"
- Rail Spacing: 2 3/8
- Designed for babies 9 months and younger, 200lb weight capacity.
- Lifetime Warranty
- Made in the USA

Options:

- Patient Weighing System
- 3" Mattress
- Crib Bumper Pads
- Under Crib Storage Bin with
- Oxygen Tank Cutouts
- IV Pole with Bracket
- Custom Powder Coat Colors

| Model # | Overall Size | Finish | Sleeping Surface |
|---------------|--------------------|--------|---------------------------|
| E1965CL-MF | 24"W x 36"L x 44"H | Epoxy | Fowler Deck |
| E1965FWCR | 24"W x 36"L x 44"H | Epoxy | Crank Fowler |
| E1965FWCR-3SD | 24"W x 36"L x 44"H | Epoxy | Crank Fowler, 3 Side Drop |
| E2065CL-MF | 24"W x 36"L x 44"H | Epoxy | ICU Series, Fowler Deck |
| E2065FWCR | 24"W x 36"L x 44"H | Epoxy | ICU Series, Crank Fowler |



Model E1965CL-MF



Model E2065FWCR