



## Economy IV Pole

*4 Hook Model: IVSTAND-04ES*



### Features:

- 4 rubber wheels provide smooth transport and maneuverability
- Removable hooks make this economy IV pole versatile for the end user
- Easily converted to 2 hook or 4 hook with easy release push-pin
- Chrome plated steel with weighted base provides strength, durability and reduces risk of tipping
- Locking collar allows for easy height adjustment
- 1 Year Warranty

### Options:



Foley Bag Hook



Also Available 2 Hook Model:  
**IVSTAND-02ES**