

## ICU Klimer Cribs

Our ICU Klimer Cribs offers the same quality construction with Vinyl side and end panels with Easy Access Safety Gate System for the ICU. Available in both chrome and epoxy finish in a variety of models to fit your facility's needs!

### Standard Features:

- Sleeping surface options: Gatch Deck, Trendelenburg Deck, Fowler Deck, Flat Pan Deck, and Gatch/Trend Deck
- Vinyl side and end panels
- Solid pan sleeping system in a non-helical design
- Easy access ICU safety gate system
- 5" Casters; 1 swivel, 1 steering and 2 braking
- Safety extended included on all gatch decks
- Available in chrome, epoxy or stainless steel finish
- Lifetime Warranty on all Steel Components  
1 Year on all other parts




# Side-by-Side Specifications

<i>ICU INFANT SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E2072CL	30"W x 44"L x 78"H	Epoxy	Flat Deck
C2072CL	30"W x 44"L x 78"H	Chrome	Flat Deck
E2072CLT	30"W x 44"L x 78"H	Epoxy	Trendelenburg Deck
C2072CLT	30"W x 44"L x 78"H	Chrome	Trendelenburg Deck
E2072FWCR	30"W x 44"L x 78"H	Epoxy	Crank Fowler Deck
C2072FWCR	30"W x 44"L x 78"H	Chrome	Crank Fowler Deck
<i>ICU CHILD SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E2082CL	30"W x 60"L x 78"H	Epoxy	Flat Deck
C2082CL	30"W x 60"L x 78"H	Chrome	Flat Deck
E2082CLT	30"W x 60"L x 78"H	Epoxy	Trendelenburg Deck
C2082CLT	30"W x 60"L x 78"H	Chrome	Trendelenburg Deck
E2082CG	30"W x 60"L x 78"H	Epoxy	Gatch Deck
C2082CG	30"W x 60"L x 78"H	Chrome	Gatch Deck
E2082CGT	30"W x 60"L x 78"H	Epoxy	Gatch/Trendelenburg Deck
C2082CGT	30"W x 60"L x 78"H	Chrome	Gatch/Trendelenburg Deck
E2082FWCR	30"W x 60"L x 78"H	Epoxy	Crank Fowler Deck
C2082FWCR	30"W x 60"L x 78"H	Chrome	Crank Fowler Deck
<i>ICU YOUTH SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E2092CG	36"W x 72"L x 78"H	Epoxy	Gatch Deck
C2092CG	36"W x 72"L x 78"H	Chrome	Gatch Deck

## Optional Accessories:

- Patient Weighing System
- Vinyl Top
- Crib Bumper Pads
- IV Pole with Bracket
- Steering Casters
- Dual Access Gates
- Mattress
- Central Locking Casters
- Undercrib Storage with O2 Cutouts
- Teething Guards
- Dialysis Gate



**NOVUM**  
Medical Products, Inc.

Phone | 716.759.7200  
Toll Free | 800.274.2742  
Fax | 716.759.0700  
Email | info@novummed.com  
www.novummed.com