

ICU Standard Cribs

Standard Features:

- Sleeping surface options: Gatch Deck, Trendelenburg Deck, Fowler Deck, Flat Pan Deck and Gatch/Trend Deck
- Safety extender included on all gatch decks
- Solid pan sleeping system in a non-helical design
- Easy access ICU safety gates
- 5 Increments of Rail Stopping Points
- 5" Casters, 1 steering, 1 swivel and 2 braking
- Rail Spacing 2 3/8"
- Crib Height: 62"
- Deck Height (from floor): 35"
- Available powder coated, chrome or stainless steel finish
- Made in the USA
- **Lifetime Warranty**

Optional Accessories:


- Patient Weighing System
- Mattress
- Cage, Bubble or Plexi Top
- Central Locking Casters
- Crib Bumper Pads
- Undercrib Storage with O2
- Cutouts
- IV Pole with Bracket
- Dialysis Gate
- Dual Access Gates



Novum ICU Hospital Cribs are designed to complement any medical setting, while ensuring the safety and comfort of the patient. We offer a selection of accessories and customizable options to meet your specific needs as well.

Side-by-Side Specifications

<i>INFANT SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E2080CL	30"W x 44"L x 61"H	Epoxy	Flat Deck
C2080CL	30"W x 44"L x 61"H	Chrome	Flat Deck
E1980CLT	30"W x 44"L x 61"H	Epoxy	Trendelenburg Deck
C2080CLT	30"W x 44"L x 61"H	Chrome	Trendelenburg Deck
E2080FWCR	30"W x 44"L x 61"H	Epoxy	Crank Fowler
C2080FWCR	30"W x 44"L x 61"H	Chrome	Crank Fowler
<i>CHILD SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E2081CL	30"W x 60"L x 61"H	Epoxy	Flat Deck
C2081CL	30"W x 60"L x 61"H	Chrome	Flat Deck
E2081CLT	30"W x 60"L x 61"H	Epoxy	Trendelenburg Deck
C2081CLT	30"W x 60"L x 61"H	Chrome	Trendelenburg Deck
E2081CG	30"W x 60"L x 61"H	Epoxy	Gatch Deck
C2081CG	30"W x 60"L x 61"H	Chrome	Gatch Deck
E2081CGT	30"W x 60"L x 61"H	Epoxy	Gatch/Trendelenburg Deck
C2081CGT	30"W x 60"L x 61"H	Chrome	Gatch/Trendelenburg Deck
E2081FWCR	30"W x 60"L x 61"H	Epoxy	Crank Fowler
C2081FWCR	30"W x 60"L x 61"H	Chrome	Crank Fowler
<i>YOUTH SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E2083CG	36"W x 72"L x 61"H	Epoxy	Gatch Deck
C2083CG	36"W x 72"L x 61"H	Chrome	Gatch Deck



NOVUM
Medical Products

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com