

IV Wall Hook

IV-WALL

IV hooks that can be easily mounted on the wall, providing ease of access that is not offered by a traditional IV pole.

Features

- Stainless Steel Hook
- Hook Diameter: 1/4"
- 16 Gauge powder coated white wall bracket
- Plastic caps are FDA compliant
- #10 screws (2) with board anchors are included
- Hook sits approximately 9.5" inches away from wall when as its longest.

Benefits

- Low Profile
- Low Cost
- Value of Product
- Ease of access compared to traditional IV pole

