

NOVUM

Medical Products

Klimer Cribs

Novum Medical Products' series of Klimer Cribs offers the same quality construction with Vinyl side and end panels. Designed to keep overactive kids safe in their cribs. Available in both chrome and epoxy finish in a variety of models to fit your facility's needs!

Standard Features:

- Sleeping surface options:
 - Gatch Deck, Trendelenburg Deck, Fowler Deck, Flat Pan Deck, and Gatch/Trend Deck
- Vinyl side and end panels
- Solid pan sleeping system in a non-helical design
- Easy access safety gates
- 5" Casters, 2 swivel and 2 braking
- Safety extended included on all gatch decks
- Available in chrome, epoxy or stainless steel finish




Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com

Side-by-Side Specifications

<i>INFANT SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E1972CL	30"W x 44"L x 78"H	Epoxy	Flat Deck
C1972CL	30"W x 44"L x 78"H	Chrome	Flat Deck
E1972CLT	30"W x 44"L x 78"H	Epoxy	Trendelenburg Deck
C1972CLT	30"W x 44"L x 78"H	Chrome	Trendelenburg Deck
E1972FWCR	30"W x 44"L x 78"H	Epoxy	Crank Fowler Deck
C1972FWCR	30"W x 44"L x 78"H	Chrome	Crank Fowler Deck
<i>CHILD SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E1982CL	30"W x 60"L x 78"H	Epoxy	Flat Deck
C1982CL	30"W x 60"L x 78"H	Chrome	Flat Deck
E1982CLT	30"W x 60"L x 78"H	Epoxy	Trendelenburg Deck
C1982CLT	30"W x 60"L x 78"H	Chrome	Trendelenburg Deck
E1982CG	30"W x 60"L x 78"H	Epoxy	Gatch Deck
C1982CG	30"W x 60"L x 78"H	Chrome	Gatch Deck
E1982CGT	30"W x 60"L x 78"H	Epoxy	Gatch/Trendelenburg Deck
C1982CGT	30"W x 60"L x 78"H	Chrome	Gatch/Trendelenburg Deck
E1982FWCR	30"W x 60"L x 78"H	Epoxy	Crank Fowler Deck
C1982FWCR	30"W x 60"L x 78"H	Chrome	Crank Fowler Deck
<i>YOUTH SERIES</i>			
Model #	Overall Size	Finish	Sleeping Surface
E1992CG	36"W x 72"L x 78"H	Epoxy	Gatch Deck
C1992CG	36"W x 72"L x 78"H	Chrome	Gatch Deck

Optional Accessories:

- Patient Weighing System
- Vinyl Top
- Crib Bumper Pads
- IV Pole with Bracket
- Steering Casters
- Dual Access Gates
- Mattress
- Central Locking Casters
- Undercrib Storage with O2 Cutouts
- Head End Safety Gate
- Dialysis Gate
- Teething Guards



NOVUM
Medical Products

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com