



Non-Tilt Overbed Table

Model NV-67D



Features

- Walnut, wood grain
- Low pressure laminate top
- Table top can be raised or lowered in infinite positions between 28"-45"
- Silver steel "H" base provides security and stability
- 15" x 30" table top
- Top locks securely when height adjustment handle is released
- Table top can be raised with slightest upward pressure
- Silver Vein Base and Mast
- 2" swivel casters allow for easy maneuverability

Specifications

- Height Adjustable: 28" - 45"
- Base Dimensions 15.5"x 26.5"
- Tabletop Dimensions: 15" x 30"
- Weight Limit: 40 lbs.