


Acute Care Beds

NV-ACB-A02



Features:

- Deck Height Range: 11" to 31.5"
- Manual CPR Quick Release Head Deck
- High Impact Quick Release Head and Foot Boards
- High Impact Side Rails
- Central Locking System with 5" Casters
- Caster Base Covers
- Drainage Hooks

Electronic:

- Pillow Speaker with Nurse Call Two Way Communication and Entertainment Controls
- Out of Bed Alarm
- Under Bed Lighting
- Battery Backup

Controls:

- Handset: 10 Button Backlit, Non-locking
 - Head, Foot, Comfort Chair, Elevation, Trendelenburg, Under Bed Light
- Y-Cable (handset plugs into either side of bed)
- Foot Board Control
 - All Handset Features Plus Battery Power Indicator

Specifications:

- Weight Capacity: 600 lbs
- Model NV-ACB-A02: 35" x 80"
- Model NV-ACB-A02-L: 35" x 84"
- Deck Height Range: 11" - 31.5"
- 5" Casters
- Self Leveling Motors
- 15 Year Warranty on Decks & Frame Welds
- 4 Year Warranty on Electronics