



Acute Care Beds



Features:

- Deck Height Range: 11" to 31.5"
- High Impact Quick Release Head and Foot Boards
- High Impact Side Rails
- Central Locking System with 5" Casters
- Caster Base Covers
- Drainage Hooks

Electronic:

- In Bed Scale with Digital Readout
- One Button CPR Quick Release Head Deck
- Pillow Speaker with Nurse Call Two Way Communication and Entertainment Controls
- Out of Bed Alarm with Nurse Call Button
- Under Bed Lighting and Backlit Handset
- Battery Backup

Controls:

- Handset: 10 Button Backlit, Non-locking
 - Head, Foot, Comfort Chair, Elevation, Trendelenburg, Under Bed Light
- Y-Cable (handset plugs into either side of bed)
- Nurse Panel: All handset features plus angle and height indicators, CPR button
- Scale Control and Readout

Specifications:

- Weight Capacity: 600 lbs
- *Model NV-ACB-A04*: 35" x 80"
- *Model NV-ACB-A04-L*: 35" x 84"
- Deck Height Range: 11" - 31.5"
- 5" Casters
- Self Leveling Motors
- 15 Year Warranty on Decks & Frame Welds
- 4 Year Warranty on Electronics