



Acute Care Electric Bed

Model NV-EAC-002



Features:

- Separate bed height and head/foot section adjustment
- 3 motors provide quiet, smooth operation
- 80" x 36" spring deck surface
- Durable powder coated epoxy finish
- Bed rests on stabilizing legs in lowest position
- Equipped with locking casters
- Sturdy, attractive head and foot boards
- 450 lb. Weight Capacity

Specifications:

- Bed Minimum Height: 15"
- Bed Maximum Height: 26"
- Overall Length: 87"
- Overall Width: 36"
- Product Weight: 179 lbs