



# Long Term Care Electric Bed

*Model # NV-ELTC-003*



## Features:

- The transition box is mounted on the foot section and allows you to change the rotation of the high-low shaft
- Wood grain panels are break and scratch resistant
- Foot section weighs 50% less than conventional foot sections to make deliveries safe and easy
- Unique motor is completely self-contained to reduce weight and noise, can be installed and removed with patient in bed and is UL approved.
- In the event of a power failure, a 9 volt battery, located in the motor, can lower the head and foot sections nine times. No crank needed.
- New and improved hand pendant has large, easy to use controls
- Channel frame construction provides superior strength and reduced weight
- Zinc coated spring deck