

Long Term Care Adult Bed



Features:

- Electric Bed: Hi/Lo, Head, Knee
Optional: Trendelenburg and Reverse Trendelenburg
- Constructed for bed to go straight up and straight down
- Designed with rib steel deck sleep surface for safety
- Bed rolls and locks at any height
- Lightweight for fast, quiet and smooth travel
- Swivel braking casters at all four corners
- Easy to clean and maintain
- Locking or Non-locking handset
- Color Options: Light Oak, Dark Oak, Wild Cherry, Walnut, Mahogany

Specifications:

- Length Options:
Model NV-ELTC-5ER76: 76"
Model NV-ELTC-5ER80: 80"
Model NV-ELTC-5ER84: 84"
- Width: 35"
- Width with side rails: 40"
- Width with assist bars: 36.5"
- Height Minimum: 7.5"
- Height Maximum: 28"
- Head Angle: 70°
- Knee Angle: 45°
- Weight Capacity: 600 lbs.
- Warranty: Electronics: 4 Years
Welding: Lifetime