




Pressure Redistribution Mattress

Model # NV-PPM7



Features:

- High quality therapeutic foam mattress provides pressure redistribution and shear/ friction reduction
- Drive's most popular therapeutic foam mattress maximizes value and comfort with 3 full layers of latex free foam
- Full width top layer uses unique die-cutting to optimize zoned pressure redistribution
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
- High density foam perimeters provide solid edge support and facilitate safe patient transfers
- Available with or without raised side rails
- Visco foam heel section offloads pressure from vulnerable heels
- Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633
- 350 lb. weight capacity

Model Numbers

NV-PPM7-80	36"W x 80"L x 6"H	Without Raised Side
NV-PPM7-80-RS	36"W x 80"L x 6"H	With Raised Side
NV-PPM7-84	36"W x 84"L x 6"H	Without Raised Side
NV-PPM7-84-RS	36"W x 84"L x 6"H	With Raised Side