



Pressure Redistribution Mattress

Model # NV-PPM8

Features:

- Outstanding combination of multi-layered die-cut and memory foam available with or without raised side rails.
- Mattress includes a 6" to 5" heel slope and a "channel cut" center layer of foam under the foot section to target heel ulcer prevention.
- Top layer uses unique die-cutting and a true memory foam heel section to optimize zoned pressure redistribution.
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress.
- Stretch Cover is fire retardant, fluid proof, tear resistant, low shear, and easy to clean.
- High Density Foam Perimeter provides solid edge support and facilitates safe patient transfers.
- Weight Capacity: 350 lbs.
- 3 Layers of Foam
- Mattress Cover Material: Stretch Fluid-Resistant Nylon Top & Waterproof Vinyl Bottom
- Fire Rating: CFR 16 Part 1633 Federal Fire Code
- Warranty: 10 Year Limited Foam, 1 Year Cover



Model Numbers

Model Number	Dimensions	Description
NV-PPM8-80	80" (L) x 36" (W) x 6" (H)	Multi-Layered/Multi-zoned Foam Mattress
NV-PPM8-84	84" (L) x 36" (W) x 6" (H)	Multi-Layered/Multi-zoned Foam Mattress
NV-PPM8-80-RS	80" (L) x 36" (W) x 6" (H)	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out
NV-PPM8-84-RS	84" (L) x 36" (W) x 6" (H)	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out

Phone | 716.759.7200
 Toll Free | 800.274.2742
 Fax | 716.759.0700
 Email | info@novummed.com
www.novummed.com