



## Eye/Oral Surgery Stretcher

*Model NV9000-EYE*



### Features:

- Four IV pole holders
- Roller bumpers, Side bumpers
- Adjustable fowler section
- 8" Tente® wheel casters, 5th wheel steering
- Fold-down siderails
- Two side mounted pedals for patient surface height adjustment
- Four-caster braking system
- Ergonomic push handles
- Head and foot-end mattress retainer
- Articulating headpiece

### Specifications:

- Overall Dimensions: 86"L X 32.5"W
- Patient Surface: 83"L X 26"W
- Weight Capacity: 550 lbs
- Height Range: 21.25" - 33"
- Head Piece:
  - Height: 4"
  - Side to Side: +/- 15°
- Positioning:
  - Head: 81.5°
  - Knee: 24°
- Trendelenburg: 15.5°
- Reverse Trendelenburg: 15.5°
- Caster Diameter: 8"
- 3 Year Limited Warranty