



Plexi Crib Top



Features:

- Allows for an unobstructed view of the crib interior
- Made of plexiglass with an aluminum frame
- Lightweight for easy lifting
- Sides slide up by squeezing two black tabs and pushing upwards to lock in place
- Can be mounted to any of Novum Medical's Standard and ICU cribs
- Halfway and full locking points when lifting up

Models/Specifications:

- **PT3044:** For 30" x 44" Infant Cribs
- **PT3060:** For 30" x 60" Child Cribs
- **PT3672:** For 36" x 72" Youth Cribs

Best used for patients that are active who may be at risk to climb out of a regular crib. The plexi crib top helps to keep the child safely in the crib while still allowing them to easily see their surroundings.