



IV Power Pole

Model IVPP-09

The Power Pole infusion stand is constructed primarily out of steel with powder coat finish. The base has heavy duty casters, was engineered to withstand 40 lbs. of fluid, and tip tested to OSHA standards.



Features

- Internally mounted gas cylinder activated by foot pedal
- Elevates 89" tall
- Heavy duty 26" base
- Removable 4 position Rams Hook
- 3 Position Adjustment Handles
- 3" Swivel Casters (6)
- Latex Free



Tired of lifting 3000ML bags?

- Elevate 4 bags with one hand
- Same size as conventional IV Tower
- No cords, plugs or batteries
- Designed with injury prevention in mind
- Perfect for arthroscopies