

Replacement Pads and Mattresses

Novum Medical Products can replace any make or model of your existing hospital mattress pad, stretcher pad, operating table pad or armboard pad with a number of intelligent options.

Foams

All foam fillers are LATEX FREE and conform to CAL 117.

Standard


High-density polyurethane foam. This pad is recommended for shorter procedures.

Softcare


This pad incorporates a layer of pressure relief Visco-Elastic foam (a.k.a. Memory Foam) to aid in the prevention of pressure ulcers during longer procedures. The balance of the foam is high-density polyurethane foam.




Standard Foam



Softcare Foam



Softcare Foam with Gel



Softcare Foam with Rubber Base

Covers

All cover materials are Latex Free, radiolucent, waterproof, stain resistant, anti-microbial, anti-bacterial, washable and conform to California Technical Bulletin 117 (CAL117).

Naugahyde


Standard conductive vinyl sheeting top and bottom cover.

Comfort

Top cover material is two-way stretch with soft, stretchable, urethane-coated, conductive material.

Fusion

Top cover material is four-way stretch with soft, stretchable, urethane-coated, conductive material.



Replacement Pads and Mattresses

Part Numbers

NV-SP-NS	Stretcher Pad, Naugahyde (vinyl) cover, Standard Foam.
NV-SP-SS	Stretcher Pad, Comfort Cover (nylon), Standard Foam.
NV-SP-FV	Stretcher Pad, Fusion Cover, Viscoe (memory) Foam.
NV-OR-NS	Operating Room Table Pad, Naugahyde (vinyl) cover, Standard Foam.
NV-OR-SS	Operating Room Table Pad, Comfort Cover (nylon), Standard Foam.
NV-OR-FV	Operating Room Table Pad, Fusion Cover, Viscoe (memory) Foam.
NV-PAD-GEL	Gel overlay for Stretcher/ OR Table pads.
NV-PAD-WELD	Welded Seam for Stretcher/ OR Table pads.
NV-PAD-RBB	Rubber Base for Stretcher/ OR Table pads.
NV-PAD-POI	Add one inch of foam for Stretcher/ OR Table pads.
NV-PAD-BFC	CAL133 (Boston Fire Code) for Stretcher/ OR Table Pads
NV-PAD-VIS	Softcare Foam for Stretcher/ OR Table pads.



Stretcher Pad



Armboard Pad



OR Table Pad



Boot Pad

Ordering Steps

Step 1)

Provide Novum with **Make and Model** of your Stretcher (SP), OR Table (OR) Pad or the dimensions.

Step 2)

Select your foam, cover and options.

Step 3)

Submit information online for a quick quote or contact Novum with information for quote.