



Restraint Straps

Features:

- Strong, cam-style metal buckles with teeth that lock securely.
- Fits tables with or without side rails.
- 2" straps are available at a variety of lengths.
- Long enough to go around the patient and table.
- Latex-free mid-panel is 18"L x 5"W.
- Available with hooks (for attaching to side-rail/ universal bar) or without hooks.
- Fits all surgical tables at a fraction of the cost of straps from table manufacturers.
- Entire strap is easy to clean: simply wipe with soap and water or disinfectant.



Item Numbers:

- NV-RS96: 96" Strap
- NV-RS108: 108" Strap
- NV-RS120: 120" Strap
- NV-RS132: 132" Strap
- NV-RSH: Restraint Strap Hook

