

iSERIES Overbed Tables

Features:



- Easy to clean with thermofoil or highly durable, 3D Laminate seamless top.
- Anodized aluminum column and heavy duty steel base with color options
- Low profile U-base for easy mobility to fit with beds, stretchers and recliners (1.75)"
- Pneumatic lift mechanism to adjust height
- Finish options available for table tops
- Optional Ultra Low Base, 1.34" to fit under beds with 5th wheels and other furnishings
- Low profile base casters



Specifications


- Overall Dimensions: 17"W x 32"L
- Height Range: 28" - 44"
- Flip Side Table: 13.5"
- Weight Capacity: 100 lbs (evenly distributed)




NOVUM
Medical Products

Finish Options

Rigid Thermoformed/3DL



Thermoformed




Premium Upgraded Finish Options—100% Seamless RTF

Request a swatch kit for additional color matching quality.



Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com

Base Color Options



Gray



Cream

Table Top Design Options



Rigid Thermo/3DL



Thermoformed



Kidney Shape



Peanut Shape

Ultra-Low Base



The lowest overbed table base in the industry!



Add-on Options




Urinal Hook



Storage Baskets



Dual Vanity with Mirror



NOVUM
Medical Products

iSERIES Casegood Line



1-Door Bedside Table

29"H x 20"W x 19"D

*Shown in...
Shaker Silhouette Style
Auburn Finish*

**Elegant.
Durable.
Affordable.**

**Build Your
Own!**

3-Drawer Bedside Table


30"H x 20"W x 19"D

*Shown in...
Shaker Style
Cocoa Pecan Finish*



iSERIES Designs
INNOVATIVE INTELLIGENT INSPIRED

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com



Thermoformed Laminate Finishes


Furniture Options

Thermoformed




Premium Upgraded Finish Options—100% Seamless RTF

Rigid Thermoformed



Textured finish/reclaimed wood look



Request a swatch kit for additional color matching quality.

Style Collections

Hardware Selections

Simplicity (1)



Shaker (2)



Shaker Silhouette (3)



Designer (A)




Archway (B)




Channel (C)



Build Your Own!



Thermoformed



Materials and Construction

- Tops and Door/Drawer Fronts: 100% Seamless RTF— $\frac{3}{4}$ " MDF composite panels with seamless rigid thermofoil laminate surface with black or white melamine backer
- Exterior Shell: $\frac{3}{4}$ " M2 particle board panels with pvc edge banding (Classic Collection also available in 100% Seamless RTF)
- Interior Shell: Fully-finished matching

Color swatches available upon request.

Drawer Boxes

- Vinyl-wrapped & moisture-resistant box
- Adjustable drawer faces
- Optional lock available for top drawer (bedside tables only)

Drawer Slides and Hinges

- Soft-close drawer slides
- Positive stop prevents tip outs
- Up to 100 lb. max weight
- Hidden European style hinges

Rigid Thermoformed



Materials and Construction

- Tops, Door/Drawer Fronts and Exterior Shell: 100% Seamless RTF— $\frac{3}{4}$ " MDF composite panels with seamless rigid thermofoil laminate surface with black or white melamine backer
- Interior Shell: Choice of black or white melamine

Tops and Door/Drawer Fronts

- Resistant to stains, moisture, and wear
- Seamless, hygienic and easy up-keep
- Recessed for spill-containment (bedside tables only)


Pulls

- ADA inspired easy-to-grip Pulls
- Choice of 3 Styles

Other Options

- Heavy-duty Casters
- Top Drawer Keyed Lock (bedside tables only)

Build Your Own!



NOVUM
Medical Products

iSERIES Casegood Line



Build Your Own!

Wardrobe

72"H x 36"W x 25"D

Shown in...
Shaker Style
Cocoa Pecan Finish

Shown in...
Simplicity Style
Ganache Finish



4-Drawer Dresser

38"H x 34"W x 19"D

3-Drawer Dresser

29"H x 32"W x 17.5"D



Headboard

20"H x 36"W

Shown in...
Shaker Silhouette Style
Cider Finish

Footboard

15"H x 36"W

iSERIES Designs
INNOVATIVE INTELLIGENT INSPIRED

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com

All recliners are made with soft, medical grade vinyl in various color options and feature spectacular lumbar support, wide shoulder room and multiple accessories to choose from.



Model IRC-T-PA10-A

Push Arm Recliner

Push on the arm to recline
3" swivel locking casters
Seat Height: 19"
Warranty: Limited 7 years

28.5"W x 35"D x 46"H
Weight Capacity: 500 lbs
Made in the USA



Model IRC-T-PA45-A

Trendelenburg Recliner

Push on the arm to recline
Foot pedal trendelenburg
3" swivel locking casters
Warranty: Limited 7 years

30.5"W x 34"D x 42"H
Seat Height: 17"
Weight Capacity: 500 lbs
Made in the USA



Model IRC-T-LO51-A

Sleeper Recliner

Lever operated footrest
Fingerpull operated back
3" swivel locking casters
Warranty: Limited 7 years

28.5"W x 35"D x 46"H
Seat Height: 19"
Weight Capacity: 500 lbs
Made in the USA



Model IRC-T-LO52-A

Sleeper Trendelenburg Recliner

Lever operated footrest
Fingerpull operated back
Foot pedal trendelenburg
3" swivel locking casters
Warranty: Limited 7 years

28"W x 34"D x 44"H
Seat Height: 17"
Weight Capacity: 500 lbs
Made in the USA



Wall Saver Recliner

Lever operated footrest
Push arm operated back
Stationary base
Warranty: Limited 7 years

28.5"W x 35"D x 46"H
Seat Height: 19"
Weight Capacity: 350 lbs
Made in the USA

Model IRC-S-LO20-A



Rocker Recliner

Lever operated footrest
Push arm operated back
Stationary base
Chair automatically stops and locks when reclined
Warranty: Limited 7 years

28.5"W x 35"D x 46"H
Seat Height: 19"
Weight Capacity: 350 lbs
Made in the USA

Model IRC-S-LO30-A



Glider Recliner

Lever operated footrest
Push arm operated back
Stationary base
Chair automatically stops and locks when reclined
Warranty: Limited 7 years

28.5"W x 35"D x 46"H
Seat Height: 19"
Weight Capacity: 350 lbs
Made in the USA

Model IRC-S-LO35-A



Push Arm Recliner

Extremely comfortable with medical grade vinyl. Available with cut-out arms, full arms and various back heights.
Seat Height: 17"
Warranty: Limited 7 years
Made in the USA

31"W x 26.5"D x 31"H
Weight Capacity: 350 lbs

Model ILC-S-VER

iSeries Color Options

Wood Arm Caps



NATURAL

CHERRY

WALNUT

Polyester Arm Caps



BLACK

GRAY

TOAST

CONCRETE

Vinyl

SNOW WHITE	SUPER WHITE	ADOBE WHITE	ALABASTER	FLINT GRAY	DOVE	PLATINUM	TOAST	PEACH	PINK	IRIS
ROSE	MAUVE	TEA ROSE	PARCHMENT	SAND	SHELL	DEWBERRY	BLUSH	BRITISH TAN	DEEP CLAY	TAUPE
GREEN TEA	DOE	SADDLETAN	OAK	TEA LEAF	CHESTNUT	MOCHA	ESPRESSO	SUN YELLOW	GOLDENROD	MANDARIN ORANGE
TOMATO	AMERICAN BEAUTY	CANDY APPLE	PAPRIKA	CLARET	BURGUNDY	RASPBERRY	SANGRIA	WINEBERRY	PLUM	WOOD VIOLET
PURPLE IRIS	GRAPE	NEW PURPLE	CONCORD	OCEAN GREY	SPACE BLUE	COLONIAL BLUE	BLUE RIDGE	MARINA	ROYAL	BLUEBERRY
REGIMENTAL BLUE	IMPERIAL BLUE	ARTICHOKE	OLIVE GREEN	VIZCAYA PALM	NORTH WOODS	IRISH SPRINGS	EMERALD	FOREST	YEW GREEN	DARK AQUA
CELADON	DUSTY JADE	TEAL	AQUA GREEN	AEGEAN	GROTTO	GUN METAL	GRAPHITE	BLACK		