


Pressure Redistribution Stretcher Pad

Model NV-PPM12



Features:

- Easy to clean, Water resistant
- Antimicrobial, antifungal, antibacterial
- Fluid Resistant and Tear Resistant
- Flammability CA 117
- Foam removes excess heat and moisture by allowing air to pass through it without losing its properties
- Memory foams combined with our high-resilience foams to create an incredibly comfortable surface
- Perfect combination for intensive care applications, where a less humid bed surface is critical.
- Pads can be custom designed to your specifications

Dimensions:

- 6" Thickness
- Available Sizes:

35"W x 76"L, 35"W x 80"L, 35"W x 84"L

39"W x 76"L, 39"W x 80"L, 39"W x 84"L

42"W x 76"L, 42"W x 80"L, 42"W x 84"L