

Medical Recliner

Model RC302



Features:

- 3 positions: upright, deep recline and elevated legrest/footrest
- Retractable lock bar secures the chair in the desired position
- Self adjusting suspension in the seat prevents bottoming out
- Side panels "pop-off" for easy cleaning and maintenance
- Large, blow molded tray locks in 5 positions and can be stored on the side of the chair
- 5" casters, 2 locking for easy maneuverability
- Wide, comfortable armrests are angled to provide support in any reclined position
- Sturdy steel frame is easy to clean and maintain
- Moisture barrier on seat
- Comfortable "built in" headrest
- Foam and vinyl meet California Flammability Codes
- Standard and Boston Fire Department available

Specifications:

Overall Height: 47"
Overall Width: 25"
Overall Depth: 41"
Seat Height: 21"
Seat Width: 19"
Seat Depth: 19"
Weight Capacity: 250 lbs

Fabric Options:



Jade



Rosewood



Blue Ridge