

Nitro DLX Rollator

HCPCS:



FIGURE A




FIGURE B



FIGURE C

- Attractive, Euro-style design
- Handle height easily adjusts with unique push button
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- 10" front casters allow excellent steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage
- Cross-brace design allows for side-to-side folding and added stability
- Durable and comfortable firm seat separates allowing unit to fold
- Handsome and removable deluxe zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight frame

PRODUCTS

Description

| | |
|---------------|--|
| RTL10266BL-HS | Nitro,FirmFoam,BlueFrame,RTL 1/ea |
| RTL10266CH-HS | Nitro,FirmFoam,ChampagneFrame,RTL 1/ea |

SPECIFICATIONS

Description

| | |
|-----------------|---|
| Handle Height | 33.5" - 38.5" |
| Seat Dimensions | 18" (W) x 10.75" (D) |
| Seat Height | 20" |
| Unit Dimensions | 28" (L) x 24" (W) |
| Warranty | Frame: Limited Lifetime Brake Handle and Housing: 5 Years |
| Weight | 22 lbs |
| Weight Capacity | 300 lbs |

ACCESSORIES

Description

| | |
|-----|---------------------------------------|
| 892 | Cane Holder, Nitro Rollator, 1/ea, CS |
|-----|---------------------------------------|