

## Viper Plus GT Wheelchair

HCPCS: K0004-1,K0004



FIGURE A



FIGURE B



FIGURE C

- Built-in seat rail extensions and extendable upholstery easily adjusts seat depth from 16" to 18" (Patent pending) (Figure A - C)
- All-Aluminum frame
- Dual Axle allows for multiple seat-to-floor height positions
- Adjustable height, flip-back, removable, padded arms
- Weighs under 35 lbs in basic configuration (excluding front rigging)
- Over-center cross brace eliminates need for seat guides Works with seating and positioning accessories
- Adjustable angle back easily adjusts from 5° - 20°
- Adjustable height back adjusts from 17" - 19" in 1" increments, adaptable for individual patient needs
- Durable, heavy-gauge flame-retardant nylon upholstery with hook and loop fastener flaps
- Quick release, 24" composite Mag-style wheels are lightweight and maintenance free
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Adjustable angle caster forks have three height adjustments and angle adjustability
- Tool-free adjustable swing-away footrests also feature easy-to-access release handles
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight
- Adjustable anti-tipper with flip-up wheels Never needs to be removed
- Comes with push-to-lock wheel locks
- Patent Pending

### PRODUCTS

### Description

PRODUCTS	Description
PLA416FBDAARAD-ELR	ViperPlsGT,16"FlpBk.AdjHtDskArm,ELR,1/ea
PLA416FBDAARAD-SF	ViperPlsGT,16"FlpBk.AdjHt.DskArm,SF,1/ea
PLA416FBFAARAD-ELR	ViperPlsGT,16"FlpBkAdjHtFullArm,ELR,1/ea
PLA416FBFAARAD-SF	ViperPlsGT,16"FlpBk.AdjHtFullArm,SF,1/ea
PLA418FBDAARAD-ELR	ViperPlsGT,18"FlpBkAdjHtDskArm,EL,R1/ea
PLA418FBDAARAD-SF	ViperPlsGT,18"FlpBk.AdjHt.DskArm,SF,1/ea
PLA418FBFAARAD-ELR	ViperPlsGT,18"FlpBkAdjHtFullArm,ELR,1/ea
PLA418FBFAARAD-SF	ViperPlsGT,18"FlpBk.AdjHtFullArm,SF,1/ea
PLA420FBDAARAD-ELR	ViperPlsGT20"FlpBkAdjHtDetDskArmELR,1/ea
PLA420FBDAARAD-SF	ViperPlsGT,20"FlpBkAdjHtDetDskArmSF,1/ea
PLA420FBFAARAD-ELR	ViperPlsGT20"FlpBkAdjHtDtFullArmELR,1/ea
PLA420FBFAARAD-SF	ViperPlsGT20"FlpBkAdjHtDetFullArmSF,1/ea
PLA422FBDAAR-ELR	ViperPlsGT22"FlpBkAdjHtDetDskArmELR,1/cs
PLA422FBDAAR-SF	ViperPlsGT,22"FlpBkAdjHtDetDskArmSF,1/cs
PLA422FBFAAR-ELR	ViperPlsGT22"FlpBkAdjHtDtFullArmELR,1/cs
PLA422FBFAAR-SF	ViperPlsGT22"FlpBkAdjHtDetFullArmSF,1/cs

SPECIFICATIONS	Description
Back Height	17" - 19"
Carton Shipping Weight	51 lbs   52 lbs   53 lbs   54 lbs
Height	36"
Length with Riggings	42"
Product Weight with 16 Seat	34 lbs
Product Weight with 20 Seat	34 lbs
Product Weight with 18 Seat	34 lbs
Product Weight with 22 Seat	35 lbs
Seat Depth	16" - 18"   18"
Seat To Floor Height	17.5" - 19.5"
Seat Width	16"   18"   20"   22"
Seat Width Closed	12.5"
Warranty	Limited Lifetime
Weight Limit with 16 , 18 and 20 Seat	300 lbs
Weight Limit with 22 Seat	350 lbs
Width Open	24"   26"   28"   30"
Arm Rest Type (s)	Flip Back, Detachable, Adj. Height Desk Arm   Flip Back, Detachable, Adj. Height Full Arm   Flip-Back, Detachable, Adj. Height Desk Arm

ACCESSORIES	Description
STDS851	Seat Belt, Velcro Close,1/ea, CS
STDS850	Seat Belt, Clasp Type,1/ea, CS
STDS855	Seat Belt,60" Auto-type,1/ea, CS
STDS856	Seat Belt,60"Velcro-type,1/ea, CS
STDS820	Telescoping IV Pole,W/C, 1/bx, CS
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, CS
STDS834	OverheadAnti-Theft Bar,Single,1/ea, CS
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, CS
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, CS
STDS807	Anti-Tips,Adj.w/Whls,Blk,1pr/bx, CS
STDS828GT	AmputeeAdaptersViperGT,1pr, CS
STDS8HK	Super Hemi Kit Viper / Viper GT, 1kit/cs, CS
STDS4056	One Arm drive for ViperGT, CS
STDS1002	Caster,8in Sentra,Viper,Astaire,1/ea, CS
STDS860	Forwrd Stablzrs/Univ.Fit , 1pr/bx, CS

PARTS	Description
STDS2369BLK	Back Uphol 22"x18",Black,1/bx, CS
STDS1A3510	Rear Wheel,VIPER GT,1/EA, CS
STDS1A3502	Frame Assist, 22" ViperPlusGT, 1/ea, CS
STDS1A3511	Quick Release Rear Axle, Viper GT, 1/ea, CS
STDS5D1012	24" Rear Wheel Bearing, C4,1 PR, CS
STDS1073	Spacers & Screws for Handrim, 8/set, CS
STDS2359BLK	Back Uphol 18"x16",Black,1/bx, CS
STDS2361BLK	Back Uphol 20"x18",Black,1/bx, CS
STDS1A3557	Backrest Cane,VIPER GT,1PR, CS
STDS2A3250	BackAngle/Height Adj.Vip.Plus,1/ea, CS
STDS2363BLK	Back Uphol 16"x18",Black,1/bx, CS
STDS1A3513	Front Axle Bolt , 1/pr, CS
STDS2A3511	Hand Rim,Viper GT,1ea, CS

## PARTS

## Description

STDS2364BLK	Seat Uphol 18"x16",Black,1/bx, CS
STDS2366BLK	Seat Uphol 20"x18",Black,1/bx, CS
STDS2368BLK	Seat Uphol 22"x18",ViperGT,1/ea, CS
STDS2356BLK-AD	Seat UpholsteryNylon,16"ViperPlusGT,1ea, CS
STDS2364BLK-AD	Seat UpholsteryNylon,18"ViperPlusGT,1ea, CS
STDS1A35DL	Adj.DeskArm,Lt.Vip.GT,1/ea, CS
STDS1A35DR	Adj.DeskArm,Rt.Vip.GT,1/ea, CS