

Cougar Wheelchair

HCPCS: K0005-1,K0005+E0951+K0195



FIGURE A



FIGURE B



FIGURE C

- Ultra-light, weighs less than 27 lbs (K0005)
- Forward folding back to reduce size
- Tool-free back height adjustable from 17"-20 ½"
- Seat depth adjustable from 16"-18"
- Multi-position back can adjust 15°
- Rigidity bar provides support
- Quick-release rear wheels (Figure A)
- Semi-Pneumatic tires on a 24" spoke wheel provide a smooth, true roll
- 3" of depth adjustment in the overall wheelbase
- Front casters are angle adjustable and height adjustable in three positions
- Arm rest can be flipped allowing 2" width adjustment
- Nylon upholstery is breathable, lightweight and comfortable
- Back folds forward for storage (Figure B - C)
- Adjustable aluminum wheel locks
- Swing-away footrests swing in and out
- Lightweight, composite foot plates with heel loops come standard
- Optional adjustable height arm (sold separately)

PRODUCTS

Description

AK516ADA-AELR	Cougar UltraLt,W/C,16"Det,T Arm,1/cs
AK516ADA-ASF	Cougar UltraLt,W/C16"Det,"T"A, 1/cs
AK518ADA-AELR	Cougar UltraLt,W/C,18"Det,"T"A,1/cs
AK518ADA-ASF	Cougar UltraLt,W/C,18"Det,"T",1/cs

SPECIFICATIONS

Description

Carton Shipping Weight	32 lbs
Height	37"
Length with Riggings	44"
Product Weight	27 lbs
Seat Depth	16"
Seat To Floor Height	17.5" - 19.5"
Warranty	Limited Lifetime
Weight Capacity	250 lbs
Width Closed	12"
Width Open	24" 26"

ACCESSORIES

Description

STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, CS
STDS601	Brake Extension,6",1/ea, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, CS

STDS831	Heel Loops, 1pr/bx, CS
STDS820	Telescoping IV Pole,W/C, 1/bx, CS
STDS824	Quick-Release Anti-Theft Rear Caster, 1/ea, CS