

Cruiser II Standard Hemi Wheelchair

HCPCS: K0002-1




FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Adjustable height, flip back, detachable desk length armrests
- 6" wheel lock extensions included
- Anti-tippers standard
- Seat belt included
- Built in seat depth extension easily adjusts from 16" 18"
- Embossed vinyl upholstery is durable, lightweight, and easy to clean
- Carbon steel frame with silver vein finish
- Composite, mag-style wheels are lightweight and maintenance free
- 8" front casters
- Push-to-lock wheel locks
- Available with swing-away footrests or elevating legrests
- Dual axle provides multiple seat to floor height positions
- Configured in hemi position out of box for a lower seat to floor height

PRODUCTS

Description

K216ADDA-ELR	Cruiser II 16" AdjHghtDskArm ElevLegRs
K216ADDA-SF	Cruiser II 16" AdjHghtDskArm Swing-Awa
K218ADDA-ELR	Cruiser II 18" AdjHghtDskArm ElevLegRs
K218ADDA-SF	Cruiser II 18" AdjHghtDskArm Swing-Awa
K220ADDA-ELR	Cruiser II 20" AdjHghtDskArm ElevLegRs
K220ADDA-SF	Cruiser II 20" AdjHghtDskArm Swing-Awa

SPECIFICATIONS

Description

Back Height	16"
Carton Shipping Weight	(16", 18", 20") 46 lbs. 50 lbs. 51 lbs.
Height	(Lowest Axle Position) 33"
Length with Riggings	42"
Product Weight without Riggings (16 , 18 , 20)	(16", 18", 20") 37 lbs. 39 lbs. 41 lbs.
Seat Depth	16" Standard, Adjustable to 18"
Seat To Floor Height	17.5" Standard, Adjustable to 19.5"
Warranty	Limited Lifetime
Weight Capacity	250 lbs
Width Closed	12"
Width Open	(16", 18", 20") 23" 25" 27"
Arm Rest Type (s)	Flip Back, Detachable, Adjustable Height Desk Arm

ACCESSORIES	Description
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, CS
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
14300	16x17.5 Gen.use adj back cush 16-21in, CS
WASL	Limb Support,Left,1/cs, CS
WASR	Limb Support, Right,1/cs, CS
STDS834	OverheadAnti-Theft Bar,Single,1/ea, CS
STDS806	Anti-Fold/Theft Device for w/c,1/ea, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, CS
STDS820	Telescoping IV Pole,W/C, 1/bx, CS
STDS821	Overhead AntiTheft w/IV Hks,1/bx, CS
STDS6008-1	Carry Bag,W/C,Oxygen,1/ea, CS
STDS1040S	Cup Holder, Swivel,1BX,6/CS, CS
STDSSHK-K3	Super Hemi Kit, Cruiser 3 & X4, 1/ea, CS