

Chrome Sport Full-Reclining Wheelchair

HCPCS: K0001/E1226-1,K0001 OR K0002+E1226+E0966+E0990
OR K0195



FIGURE A

- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Arms have support extension for reclining position
- Comfortable, lightweight, easy-to-clean vinyl upholstery
- Attractive chrome frame
- Dual axle position for the rear Mag-Style wheels
- Lightweight, maintenance-free, Composite Mag-Style wheels
- 8" front casters
- Padded armrests provide added patient comfort
- Wheels are positioned at the back of the frame for additional stability
- Added rigidity bar keeps the back from hammocking
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Rear anti-tippers standard
- Headrest extension with cushioned head immobilizer standard
- Comes standard with swing-away elevating leg rests
- Metal foot plates
- Comes with push-to-lock wheel locks

PRODUCTS

Description

CS16RBDDA	Chrome Sport 16" Recl Back Desk Arms 1ea
-----------	--

SPECIFICATIONS

Description

Back Height	Without Headrest: 22" With Headrest: 32"
Carton Shipping Weight with 16" Seat	65.45 lbs
Carton Shipping Weight with 18" Seat	66.45 lbs
Carton Shipping Weight with 20" Seat	67 lbs
Height	51"
Length with Riggings	49"
Product Weight with 16" Seat	56.63 lbs
Product Weight with 20" Seat	58.63 lbs
Product Weight with 18" Seat	57.63 lbs
Seat Depth	17"
Seat To Floor Height	18" / 20"
Seat Width	16" 18" 20"
Warranty	Limited Lifetime
Weight Limit with 16" Seat	300 lbs
Weight Limit with 18" Seat	300 lbs

SPECIFICATIONS	Description
Weight Limit with 20 Seat	350 lbs
Width Closed	12"
Width Open	24" 26" 28"
ACCESSORIES	Description
STDS4Y4712	Anti Tippers for SSP18RBDDA,1/ea, CS
STDS831	Heel Loops, 1pr/bx, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, CS
STDS806	Anti-Fold/Theft Device for w/c,1/ea, CS
STDS851	Seat Belt, Velcro Close,1/ea, CS
STDS856	Seat Belt,60"Velcro-type,1/ea, CS
STDS855	Seat Belt,60"" Auto-type,1/ea, CS
STDS850	Seat Belt, Clasp Type,1/ea, CS
PARTS	Description
STDS3J2424	Wheel24in Blk.Handrl,Cruiser3,SSP,BS,1/e, CS
STDS3J2413	Rear Axle Bearing,Cruiser3,1/ea, CS
STDS3J2410R	Footplate&HeelLoop,Right,Cruiser,SSP,CS, CS
STDS8M4601	ForkStemBearing,SSP2,Cruiser,1/ea, CS
STDS3J4511	BackUphol.Screw&Washer 1/ea, CS
STDS3J3900	Front Caster,8"TR39,1/ea, CS
STDELR-TF	Swingway Elev.Leg Rest w/ToolFree,1pr/cs, CS
TRHG	Hand Grips/Universal,1pr/bx, CS
STDS3006	Cable,Back Release,Reclining w/c, 1/ea, CS
STDS004609-9	Hand Rim Screw Set of 9,1Bag, CS
STDS4S4809	SpreaderBar,16",CS,Senra,Viper Rec,1/ea, CS
STDS4S4810	SpreaderBar,18",CS,Senra,ViperRec,1/ea, CS
STDS4S4811	SpreaderBar,20",CS,Senra,Viper,Rec,1/ea, CS
STDS4S4829	Anti-Tippers, CS Recliner, 1/pr, CS
STDS4S2432	Caster with 11mm Bearings, 1/ea, CS
STDS4Y4702L	ReleaseHndle w/Cable,SSP18RBDDA,LFT,1/ea, CS
STDS4Y4715	Anti Folding ADJ Knob,SSP18RBDDA,1/ea, CS
STDS8M4631	FrontAxleBolt,SS2,K3,CX4,1/ea, CS
STDS4S4823	Seat Upholstery, 16", CS Recliner, 1/e, CS
STDS4S4817R	Wheel Lock, Right, CSRB,SSPRB,1/ea, CS