

Pediatric Viper Plus Reclining Wheelchair

HCPCS: E1236 E1226-1, E1236, E1226



FIGURE A




FIGURE B

- State-of-the art hydraulic reclining mechanism
- Backrest reclines from 90° to 180°
- Double-coated chrome carbon steel frame for maximum durability
- Durable flame retardant nylon upholstery
- Padded flip-back detachable desk armrests
- Composite lightweight 24" Mag-style rear wheels with urethane tires and 6" front rubber tires
- Offset rear wheel position provides additional stability and evenly distributes users weight in all reclining positions
- Precision sealed front and rear wheel bearings
- Front caster forks are adjustable in three positions
- Standard dual axle provides easy transition of seat to hemi-level
- Push-to-lock wheel locks
- Spreader bar functions as a back stabilizer and a push handle for caregiver
- Standard with height adjustable swing-away elevating leg rests
- Durable and lightweight composite footplates
- Flip-up anti-tippers with wheels

PRODUCTS

Description

PL412RBDDA ViperPlusRecl,12"DetDskArm,1/cs

PL414RBDDA Viper,14"Recl.FlpBkDsk1/cs

SPECIFICATIONS

Description

Back Height 26"

Overall Height 40.5"

Overall Length w/ Riggings 38"

Seat To Floor Height 15", 16"

Warranty Limited Lifetime

Weight Capacity 250 lbs

ACCESSORIES

Description

STDS850 Seat Belt, Clasp Type,1/ea, CS

STDS851 Seat Belt, Velcro Close,1/ea, CS