

Bariatric Sentra EC Heavy-Duty, Extra-Extra-Wide Wheelchair

HCPCS: K0007-1,K0007



- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable reinforced nylon upholstery
- Steel, spoke wheels are strong and maintenance free
- Heavy-duty, flat-free 8" front casters are adjustable in two positions
- Extra-heavy-duty front forks
- Reinforced side frame, caster journals and dual cross bars
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Extra-large, padded armrests provide comfort and support
- Riggings have extra-large, heavy-duty, aluminum foot plates that provide maximum support, are attractive and crack proof (Must be ordered separately)
- Comes with push-to-lock wheel locks
- Weight Capacity 700 lbs

PRODUCTS	Description
STD26ECDDA	Sentra EC,26"Det.DeskArm,1/cs
STD26ECDF A	Sentra EC,26in Det.FullArm,1/cs
STD28ECDDA	Sentra EC,28"Det.DeskArm,1/cs
STD28ECDF A	Sentra EC,28"Det.FullArm,1/cs
STD30ECDDA	Sentra EC,30"Det.DeskArm,1/cs
STD30ECDF A	Sentra EC,30"Det.FullArm,1/cs

SPECIFICATIONS	Description
Back Height	18"
Carton Shipping Weight	68 lbs 70 lbs 72 lbs
Height	35"
Length with Riggings	43"
Product Weight with 26 Seat	62 lbs
Product Weight with 28 Seat	64 lbs
Product Weight with 30 Seat	66 lbs
Seat Depth	20"
Seat To Floor Height	19.5"
Warranty	Limited Lifetime
Weight Capacity	700 lbs
Width Closed	14" 14.5" 15"
Width Open	34.5" 37" 39"

ACCESSORIES**Description**

STDS829	Anti Tippers For 26",28"30"W/C,1pr/bx, CS
STDS801	Brake Handle Ext.,8", 1/ea, CS
PH-SF	Footrests, SentraHD EC XW, 1/pr, CS
STDS856	Seat Belt,60"Velcro-type,1/ea, CS
PHELR	ElevLegRst,Sentra EC HD XW,1/pr, CS
STDS820	Telescoping IV Pole,W/C, 1/bx, CS
STDS803	Oxygen Tank Holdr for W/C,1/cs, CS
STDS855	Seat Belt,60" Auto-type,1/ea, CS
STDS831	Heel Loops, 1pr/bx, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, CS

PARTS**Description**

STDSFLVB	UniversalWheelChairVinylFullArmPad,1/ea, CS
STDS1M222R	DeskArmAssmRight.,SentraECXW,1/ea, CS
STDS1M222L	DeskArmAssm.Left,SentraECXW,1/ea, CS
STDS1M223R	FullArmAssm.Right,SentraECXW,1/ea, CS
STDSY22SFR	Side Frame,Right,Sentra ECXW,1ea, CS
STDS1M223L	FullArmAssm.Left,SentraECXW,1/ea, CS
STDSY22SFL	Side Frame,Left,Sentra ECXW,1ea, CS
STDS1M2206	Back Upholstery, 26",Sentra ECXW, 1/ea, CS
STDS1M2208	BackUpholstery 28",SentraECXW,1/ea, CS
STDS1M2213	BackUpholstery30",SentraECXW,1/ea, CS
STDS1M2201	Handgrip,SentraECXW,1/ea, CS
STDS1M2200	24"SpokeWheel Assy. Complete,EC XW,1/ea, CS
STDS1M2203	RearAxleBolt,SentraECXW,1/ea, CS
STDS1M22X6	CrossFrame,28",SentraECXW,1/ea, CS
STDS1M2212	Front Caster, Sentra ECXW,1/ea, CS
STDS1M2210	FrontAxleBolt,SentraECXW,1/ea, CS
STDS1M2215	ForkStemBearingAssm.,SentraECXW,1/ea, CS
STDS1M2214	FrontFork,SentraECXW,1/ea, CS
STDS1M2226	Seat Upholstery, 26" Sentra ECXW,1/ea, CS
STDS1M2228	SeatUpholstery28",SentraECXW,1/ea, CS