

Titan P22

HCPCS: K0822-1, K0822



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Height-adjustable swivel seat features fold-down backrest
- Includes interchangeable color panels in red and blue
- Seat back reclines and features adjustable headrest
- Adjustable-length controller mount (Figure B)
- Armrests are padded, adjustable and removable
- Available with 18" and 20" wide Pan Seats (Figure A)
- Optional swing-away controller arm (Item #AA4100) sold separately
- Solid, flip-up footplate with adjustable height, depth and angle
- Positioning belt included
- Optional elevating leg rests (requires Item #AE2500 & Item #LK3JELR) sold separately
- Programmable controller
- Flat-free, non-marking tires (Figure C)
- Easy free wheel operation (Figure D)
- Anti-tip wheels

PRODUCTS

Description

TITAN18P22	Titan P22 K0822 FWD 18" x 18" Pan Seat
TITAN20P22	Titan P22 K0822 FWD 20" x 18" Pan Seat

SPECIFICATIONS

Description

Base Weight	68 lbs
Batteries (pair)	12V x 35AH (U1)
Battery Weight (pair)	44 lbs
Brakes	Electromagnetic
Caster Wheels	6" x 2" Solid
Charger	3A Offboard
Climbing Angle	6°
Controller	Dynamic Shark 40A
Dimensions	40" (L) x 24" (W)
Drive Wheels	10" x 3" Flat-free
Floor to Seat Height	20.5" - 23"
Ground Clearance	2.5"
HCPCS	K0822
Maximum Range	15 miles
Motors (2)	24V x 200W x 3900 rpm

SPECIFICATIONS	Description
Seat Weight	43 lbs
Top Speed	4 mph
Total Weight (excl. hand control)	158 lbs
Turning Radius	25"
Warranty on Batteries	6 months
Warranty on Drive Train Components	24 Months
Warranty on Electronic Controller	24 Months
Warranty on Frame	Limited Lifetime
Weight Capacity	300 lbs
ACCESSORIES	Description
AE2500	ELR Bracket - Hemi Spacing 1/EA, CS
AA4100	SwingAwayControllerArm,Titan, 1/ea, CS
AB1000	Armrest Bag, CS
AZ0060	Cup Holder, CS
PARTS	Description
TITAN-14	FootrestDepthAdjustmentBolt, Titan, 1/e, CS
TITAN-29	COMPLETESEATASSMBLY,16X16X18,TITAN,1/ea, CS
TITAN-32	COMPLETESEATASSMBLY,20X20X18,TITAN,1/ea, CS
DK-REMD01	Overlay, Dynamic LPW182001 Joystick, 1ea, CS
TITAN-31	COMPLETESEATASSMBLY,18X16X18,TITAN,1/ea, CS
TITAN-33	COMPLETESEATASSMBLY,22X19X19.5,TITAN,1/e, CS
TITAN-38	Seat Post, Titan, 1/ea, CS
LPW182001-03	ShroudtoSharkCntlrcble,CobaltW/SharkJoys, CS
C07-072-01001	Seat Plate, Titan, VW Ser. #, 1/ea, CS
VENT-86	Shroud Knob, Ventura, 1/ea, CS
LK3JELR	ELR Set,Black,1Pair eq EA, CS
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, CS
LRM412202-HC	Hub Cap, Spitfire, 1/ea, CS
S31052	Armrest Knob POV 1/ea, CS
TITAN-01	FootplateW/AttachingBracket,Titan,1/ea, CS
LDR410011-12	Armrest Knob - .5", CS
TITAN-03	Caster W/Bearings, Titan, 1/ea, CS
TITAN-05	Frame, Titan, 1/ea, CS
TITAN-04	Caster Axle, Titan, 1/ea, CS
TITAN-07	Battery Tray , Titan, 1/ea, CS