

Titan X16

HCPCS: K0816-1, K0816



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Height-adjustable swivel seat features fold-down backrest (Figure A)
- Includes interchangeable color panels in red and blue
- Seat back reclines and features adjustable headrest
- Adjustable-length controller mount (Figure B)
- Armrests are padded, adjustable and removable
- Available with 16", 18", 20" and 22" wide Captains Seat
- Optional swing-away controller arm (Item #AA4100) sold separately
- Solid, flip-up footplate with adjustable height, depth and angle
- Positioning belt included
- Optional elevating leg rests (requires Item #AE2500 & Item #LK3JELR) sold separately
- Programmable controller
- Flat-free, non-marking tires (Figure C)
- Easy free wheel operation (Figure D)
- Anti-tip wheels

PRODUCTS

Description

TITAN1616X16	Titan X16 K0816 FWD 16" x 16" Capt.Seat
TITAN1618X16	Titan X16 K0816 FWD 16" x 18" Capt.Seat
TITAN1816X16	Titan X16 K0816 FWD 18" x 16" Capt.Seat
TITAN18CSX16	Titan X16 K0816 FWD 18" x 18" Capt.Seat
TITAN2020X16	Titan X16 K0816 FWD 20" x 20" Capt.Seat
TITAN20CSX16	Titan X16 K0816 FWD 20" x 18" Capt.Seat
TITAN22CSX16	Titan X16 K0816 FWD 22" x 20" Capt.Seat
TITANLB18CSX16	Titan X16 K0816 FWD 18" x 18" Capt.Seat

SPECIFICATIONS

Description

Base Weight	68 lbs
Batteries (pair)	12V x 35AH (U1)
Battery Weight (pair)	44 lbs
Brakes	Electromagnetic
Caster Wheels	6" x 2" Solid
Charger	3A Offboard
Climbing Angle	6°
Controller	Dynamic Shark 40A
Dimensions	40" (L) x 24" (W)
Drive Wheels	10" x 3" Flat-free

SPECIFICATIONS	Description
Floor to Seat Height	20.5" - 23"
Ground Clearance	2.5"
HCPCS	K0816
Heaviest Piece Disassembled	55 lbs
Maximum Range	15 miles
Motors (2)	24V x 200W x 3900 rpm
Seat Weight	43 lbs
Top Speed	4 mph
Total Weight (excl. hand control)	158 lbs
Turning Radius	25"
Warranty on Batteries	12 Months
Warranty on Electronic Controller/Drive Train Components	14 months
Warranty on Frame	Limited Lifetime
Weight Capacity	300 lbs
ACCESSORIES	Description
AB1000	Armrest Bag, CS
AZ0060	Cup Holder, CS
AA4100	SwingAwayControllerArm,Titan, 1/ea, CS
PARTS	Description
LRC402414	Seat Bracket - Captain Seat, CS
ST-HEADREST	Headrest - Captain Seat Assembly, CS
LDM400013	Caster Fork,Titan,Intrepid,1/ea, CS
LRC402414-13	Seat Post - Captain Seat, CS
LDM400001-07	Controller Tube Knob, CS
LDR410110	Armrest Pad - Titan,Intrepid,1/ea, CS
LDR410011L	Left ArmrestAssembly,Titan,Intrepid,1/ea, CS
LDR410011R	RightArmrestAssembly,Titan,Intrepid,1/ea, CS
LRA402209	3.5 Ah Battery Charger - POV, CS
LPW182001	SharkJoystick,WCat450,Med450,Titan,1/ea, CS
LPW182001-02	SharkJoystickToShroudCable,LPW182001,1ea, CS
LDR410017	Seatbelt, CS
LDR4110	DustcapTitan,Med.450, 1/ea, CS
LDR420009-01	Anti-tipWheel-CoblT,Osprey,TitanS#1L,1/e, CS
LRA402307	50 Amp Circuit Breaker, CS
LRA402208	35AH (U1) Battery 1/ea, CS
LDM400001-04	Controller Tube - Square, CS
LDM400001-01	JoystickMountingScrewsACMPower,1pr, CS
TITAN-27	Repl.Tip for Clutch Release Hndle, Titan, CS
TITAN-28	Nut & Washer for Drive Wheel, Titan, 1ea, CS