

Nitro Aluminum Rollator, 10" Casters

HCPCS:



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E




FIGURE F



FIGURE G



FIGURE H

- Adjustable contoured backrest
- Cross-brace design for added safety
- Clip hold and side-by-side folding
- Large and convenient storage bag included
- Weight capacity: 300 lbs.
- Weighs only 17.5 lbs., easy-to-transport
-
-

PRODUCTS	Description
RTL10266	Nitro Rollator, Red Frame, 1/ea
RTL10266BK	Nitro Rollator, Black frame,1EA
RTL10266WT	Nitro Rollator, White frame,1EA

SPECIFICATIONS	Description
5 Year Limited Warranty	(brake and handle housing)
Carton Shipping Weight	21.5 lbs.
Handle Height	33.5" - 38.25"
Seat Dimensions	18" (W) x 20.5" (H) x 10" (D)
Unit Dimensions	23" (W) x 27.75" (L)
Warranty Limited Lifetime	(frame)
Weight	17.5 lbs.
Weight Capacity	300 lbs

ACCESSORIES	Description
10266-CH	Nitro Cup Holder, 1/ea, EA

PARTS	Description
1026602	HandGripLeft for Nitro,1/ea, CS
1026601-RN	HandBrakeAssembly,Left,10266,1ea, CS
1026603	HandGripRight for Nitro,1/ea, CS
1026601-W	HandBrakeAssembly,Left,10266WT,1/ea, CS
1026601-G	HandBrakeAssembly,Left,10266BK/10266BK-T, CS
1026630N	HandlebarForNitro/NitroTall,Side Mount, CS
1026604-B	BackrestAssemblyfor 10266 & 10266WT,1/ea, CS
1026605-W	BackrestLockingKnobW/Boltfor10266WT,1/pr, CS
1026605-G	BackrestLockingKnobW/Boltfor10266BK,1/pr, CS

PARTS

Description

1026604-G	Backrest Assembly for 10266BK,1/ea, CS
1026605-R	BackrestLockingKnobW/Boltfor10266,1/pr, CS
1026624	RightForkCompleteW/CasterforNitro, 1/ea, CS
1026625	LeftForkCompleteW/CasterforNitro, 1/ea, CS
1026617	FrontWheelWithBearingsforNitro, 1/ea, CS
1026616	WheelBrakeRightW/Hdwr,Nitro,Metal,1/ea, CS
1026615	WheelBrakeLeftW/Hdwr,Nitro,metal,1/ea, CS
1026612-R	RearWheelWithBearingsforNitro,Right,1/ea, CS
1026612-L	RearWheelWithBearingsforNitro,Left,1/ea, CS
1026609-FSW	SeatLiftingStrapW/AttachHrdwrefor10266WT, CS
1026609-FSG	SeatLiftingStrapW/AttachHrdwrefor10266BK, CS