

Scout 3

HCPCS: K0800-1, K0800



FIGURE A




FIGURE B



FIGURE C




FIGURE D




FIGURE E



FIGURE F

- Quick and easy disassembly
- Interchangeable color panels hide scuffs and scratches
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy free wheel operation
- Available with 12AH or 20AH batteries
- 20AH batteries increase scooter's maximum range

PRODUCTS

Description

SFSCOUT3	Scout 3-Wheel Scooter, Red/Blue, 1/ea
SFSCOUT3-EXT	Scout 3-Wheel 20AH Battery, Red/Blue

SPECIFICATIONS

Description

Base Weight	49 lbs
Batteries (pair)	12V x 12AH, 12V x 20AH (EXT Model Only)
Battery Weight (pair)	20 lbs
Brakes	Electromagnetic
Charger	2A Offboard
Climbing Angle	6°
Controller	CTE 45A
Dimensions	42.5" (L) x 20.5" (W)
Floor to Seat Height	17" - 19"
Front Wheels	8" x 2" Flat-free
Ground Clearance	2.5"
HCPCS	K0800
Heaviest Piece Weight	29 lbs
Maximum Range	9 miles (12AH) 15 miles (20AH)
Motor	24V x 270W x 4700 rpm
Rear Wheels	8" x 2" Flat-free

SPECIFICATIONS	Description
Seat Dimensions	16.5" (W) x 16" (D)
Seat Weight	16 lbs
Top Speed	4.25 mph
Total Weight (excl. hand control)	85 lbs
Turning Radius	45.5"
Warranty on Batteries	12 Months
Warranty on Drive Train Components	24 Months
Warranty on Electronic Controller	24 Months
Warranty on Frame	Lifetime
Weight Capacity	300 lbs
Wheel Base Dimensions (L x W)	32.7" x 16.5"

ACCESSORIES	Description
SF8010	Oxygen Holder Set, L Type, 1/ea, CS
SF8030	Crutch Holder Set, SType, 1/ea, CS
SF8020	Rear Basket Set, LType, 1/ea, CS
SF8050	Sun Shade for Scooter, 1/ea, CS

PARTS	Description
SC31053	Transaxle W/ Motor & Brake, Scout, 1/ea, CS
SC31052	Controller, 45A, Scout 3-Wheel, 1/ea, CS
SC31042-D	Blue Decorations, Scout 3 Deluxe, 1/ea, CS
SC31040	Front Frame, Scout 3, 1/ea, CS
SC31041-D	Red Decorations, Scout 3 Deluxe, 1/ea, CS
SC31042	Blue Decorations, Scout 3, 1/ea, CS
SC31039	Front Shroud, Scout 3, 1/ea, CS
SC31041	Red Decorations, Scout 3, 1/ea, CS
SC31044-D	Front Wheel, Scout 3 Deluxe, 1/ea, CS
SC31046	Front Axle, Scout 3, 1/ea, CS
SC31038	Front Fork, Scout 3, 1/ea, CS
SC31043	Front Wheel, Black, Scout 3, 1/ea, CS
SC31047	Battery Box Hold-Down, Scout 3 & 4, 1/e, CS
SC31049	Rear Axle Bolt Assmby, Scout 3 & 4, 1e, CS
SC31033	ConnectionBlockSet,Rear&Front,Scout,1/ea, CS
SC31026	WireHarness,Scout,CurtisElectronics,1/ea, CS
SC31012	Rear Frame, Scout, 1/ea, CS
SC31011	Rear Shroud, Scout, 1/ea, CS
SC31025	Hub Cap, Scout, 1/ea, CS
SC31037	Electromagnetic Brake, Scout, 1/ea, CS