

Spitfire DST 4-Wheel



HCPCS:




FIGURE A




FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Dynamic Suspension Technology (DST) features front and rear suspension coils for the most comfortable ride on all terrain (Figure E)
- Charge port located on the control panel and battery pack for easy charging
- Headlight
- Interchangeable color panels in red and blue
- Quick & easy disassembly
- Armrests are padded, contoured, width and angle adjustable (Figure G)
- Delta tiller with ergonomic throttle control (Figure F)
- Armrest covers with reflectors (Figure G)
- Large plastic basket with carry handle
- Height adjustable swivel seat with fold down backrest
- Anti-tip wheels

PRODUCTS

Description

SFDST4	Spitfire DST 4-Wheel Travel Scooter
--------	-------------------------------------

SPECIFICATIONS

Description

Anti-Tip Wheels	2" x 2"
Back Height	14"
Base Weight	70 lbs.
Batteries (pair)	20AH
Battery Weight (pair)	30 lbs. (20 AH)
Brakes	Electromagnetic
Charger	2A
Controller	CURTIS 45A
Dimensions	42.5" (L) x 21.25" (W) x 36" (H)
Estimated Range	15 miles
Freewheel Mode	Yes
Front Wheels	8" x 2"
Ground Clearance	2.5"
Heaviest Piece Weight	Front Section: 43 lbs.
Max Climbing Angle	6°
Motor	270W x 4700RPM
Rear Wheels	8" x 2"
Seat Dimensions	17.5" (W) x 16" (D)
Seat Weight	23 lbs.

SPECIFICATIONS	Description
Top Speed	4.25 mph
Total Weight (excl. hand control)	123 lbs.
Turning Radius	55"
Warranty on Batteries	12 Months
Warranty on Frame	Limited Lifetime
Warranty on Electronics, Controller, Drive-Train Components	24 months
Weight Capacity	325 lbs
Weight Capacity (13046 13046-HD)	325 lbs.

PARTS	Description
C11-076-00100	BatteryIndicator,LED,Scout/PhoenixHD,1ea, CS
C08-069-06100	Front Shroud, ScoutDST4, 1/ea, CS
C01-069-01701	Front Frame, ScoutDST4, 1/ea, CS
SC31052-4	Controller, 45A, Scout 4-Wheel, 1/ea, CS
DST-004R	ShroudDecorationSet,Red,ScoutDST4, 1/ea, CS
DST-004B	ShroudDecorationSet,Blue,ScoutDST4,1/ea, CS
C03-075-00600	L Steering Axle, ScoutDST4, 1/ea, CS
DST-011	Tie Rod Assembly, ScoutDST4, 1/ea, CS
C03-075-00301	Steering Knuckle, Left, ScoutDST4, 1/ea, CS
C03-075-00401	Steering Knuckle, Right, ScoutDST4, 1/e, CS
C03-075-00703	Front Fork, ScoutDST4, 1/ea, CS
C05-060-00300	Front Wheel Assembly, ScoutDST4, 1/ea, CS
C08-069-05800	Rear Shroud Handle, ScoutDST, 1/ea, CS
DST-012	Frame Latch Assembly, ScoutDST, 1/ea, CS
C04-075-00200	Rear Spring, ScoutDST, 1/ea, CS
C01-069-01901	Rear Frame, ScoutDST, 1/ea, CS
C03-075-01000	Front Spring, ScoutDST, 1/ea, CS
DST-010	ForkStmBearingSet,ScoutDST/TITANLTE,1/pk, CS
DST-008	Wire Harness, Complete, ScoutDST, 1/ea, CS
C13-069-00100	Headlight Assembly, ScoutDST, ,1/ea, CS