

Spitfire EX2 3-Wheel Travel Scooter with 12AH Battery

HCPCS: K0800



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F




FIGURE G




FIGURE H

- 6 sets of interchangeable color panels in red, blue, silver, black, white and pink
- Swivel mount cup holder included
- Extra bright LED headlight
- Charge ports located on the tiller and battery pack for easy charging
- Height adjustable swivel seat
- Adjustable padded armrests with reflectors
- Large plastic basket with removable liner
- Disassembles into 5 easy to carry pieces
- Delta tiller for easy throttle control

PRODUCTS

Description

SFEX2317FS-12	Spitfire EX2, 3-Wheel 17" Seat, 12AH
---------------	--------------------------------------

SPECIFICATIONS

Description

Back Height	16"
Base Weight	60 lbs.
Batteries (pair)	12AH 21AH
Battery Charger	2A Offboard
Battery Weight (pair)	22 lbs. (12AH) 28 lbs. (21AH)
Brakes	Electromagnetic
Controller	45A
Freewheel Mode	Yes
Front Wheels	7.5" x 2"
Ground Clearance	3"
Max Climbing Angle	6°
Max Speed	5 mph
Maximum Range	9 Miles 15 Miles
Motor	24V x 270W x 4700 rpm
Overall Dimensions	42" x 20" x 36"
Rear Wheels	8.5" x 3"
Seat Depth	16"
Seat To Floor Height	19" - 21"
Seat Weight	17 lbs.
Seat Width	17"
Total Weight (excl. hand control)	99 lbs. (12AH) 105 lbs. (21AH)
Turning Radius	46"

SPECIFICATIONS**Description**

Warranty on Batteries	12 Months
Warranty on Electronics	24 Months
Warranty on Frame	Limited Lifetime
Weight Capacity	300 lbs
Weight Capacity (13046 13046-HD)	300 lbs.

ACCESSORIES**Description**

SF8030	Crutch Holder Set, SType, 1/ea, CS
SF8010	Oxygen Holder Set, L Type,1/ea, CS
SF8020	Rear Basket Set, LType, 1/ea, CS