

## Spitfire DST Travel Scooter



HCPCS:




FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E

- Dynamic Suspension Technology (DST) features front and rear suspension coils for the most comfortable ride on all terrain
- Charge port located on the control panel and battery pack for easy charging
- Headlight
- Interchangeable color panels in red and blue
- Quick & easy disassembly
- Armrests are padded, contoured, width and angle adjustable.
- Delta tiller with ergonomic throttle control
- Armrest covers with reflectors
- Large plastic basket with carry handle
- Height adjustable swivel seat with fold down backrest
- Anti-tip wheels

### PRODUCTS

### Description

SFDST3	Spitfire DST 3-Wheel Travel Scooter
--------	-------------------------------------

### SPECIFICATIONS

### Description

Anti-Tip Wheels	2"
Base Weight	60 lbs.
Batteries (pair)	20AH
Battery Charger	2A
Brakes	Electromagnetic
Controller	CURTIS 45A
Estimated Range	15 miles
Freewheel Mode	YES
Front Wheels	8" x 2"
Max Climbing Angle	6°
Motor	270W x 4700RPM
Overall Dimensions	42.5" X 21.25" x 36" (L x W x D)
Rear Wheels	8" x 2"
Seat Dimensions	17.5" X 16" X 14" (L x W x D)
Top Speed	4.25 mph
Turning Radius	45"
Weight Capacity	325 lbs
Weight Capacity (13046   13046-HD)	325 lbs.