

Go-Lite Bariatric Steel Rollator

HCPCS: E0149 / E0156



FIGURE A



FIGURE B

- 8" casters with soft-grip tires are ideal for indoor and outdoor use
- Soft padded oversized seat
- Comes with large 20.5" (W) x 4.5" (H) x 9" (D) basket that can be mounted under seat
- Tool-free removable padded backrest for comfort
- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure A)
- Handles easily adjustable in height with self-threading knob (Figure B)

PRODUCTS	Description
10215BL-1	Rollator,Overszd w/LpBk BL,1/cs
10215RD-1	Rollator,Ovrszd,w/LpBk Red1/cs

SPECIFICATIONS	Description
Carton Shipping Weight	30.5 lbs.
Handle Height	35.25" - 39.5"
Seat Dimensions	18" (W) x 22" (H) x 13" (D)
Unit Dimensions	30.5" (W) x 26.75" (L)
Warranty	Limited Lifetime
Weight	26 lbs.
Weight Capacity	500 lbs

PARTS	Description
10215RW	RearWheel for 10215,1/ea, CS
10215F	Front Fork For 10215,1/ea, CS
10215FA	Front Axle Bolt,10215,1/ea, CS
9502P21503	Frame Fork Fastener,10215,1/ea, CS
10215W	Wheel for 10215 Rollator & C4 W/C,1/ea", CS
9501021501B	Folding Hinge Assm, BLACK,10215,1/ea, CS
950000217-ER	Wheel Brake,GoLite,Rt,1/ea, CS
10215RA	Rear Axle Bolt,10215,1/ea, CS
950000217-DL	Wheel Brake for GoLite,Lft,1/e, CS
STDS002050	Rear Axle Bearing for TR & FW, 1/ea, CS
9500021507	Handle Bars for 10215 Rollator, 1/pr, CS
10208GRP	Foam Grip,10208 Rollator,1/ea, CS
95010215JHB	HandleBar,10215J,1pair, CS
10208C	CableClip for RollatorBrake,1/ea, CS
9500004101	Cable,Fastner,1/ea, CS

PARTS	Description
9501A75000	Back Release Button,1/ea, CS
10208K	Knob,HandleHt.for Rollator,LoopBrk,1/ea, CS
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, CS
10215BSK	Basket for 10215 Rollator,1/ea, CS
10208BRAKE	Brake w/Cable,10208,1/ea, CS