

Deluxe, Trigger Release Folding Walker

HCPCS: E0135



FIGURE A



- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity (Figure A) Trigger release can be reached with hands of all sizes
- Durable, Composite Trigger Release resists cracking
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- Patented

PRODUCTS

Description

10224-4	Walker,Dual,TriggerRel,Adlt,4/cs
10225-4	Walker,DualTriggerRel.Jun.,4/cs

SPECIFICATIONS

Description

Carton Shipping Weight	27 26 lbs.
Overall Width	23.5"
Warranty	Limited Lifetime
Weight	6.5 lbs. 6 lbs.
Weight Capacity	350 lbs
Width	18" inside hand grip 20.75" inside back legs