

Folding Junior Walker, Two Button with 5" Wheels

HCPCS: E0143-1,E0143



- Easy push-button release mechanisms
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter construction ensures maximum strength
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces
- Vinyl-contoured hand grip
- Comes with wheels
- Limited Lifetime Warranty

PRODUCTS

Description

10247-4	Walker,Steel, Junior 2Button 5" Whs,4/CS
---------	--

SPECIFICATIONS

Description

Carton Shipping Weight	30.64 lbs
Depth (opened at base)	18"
Height	25" - 32"
Overall Width	24"
Warranty	Limited Lifetime
Weight	8 lbs
Weight Capacity	350 lbs
Width	17" inside hand grip 20" inside back legs