

Bariatric Aluminum Folding Walker- Two Button

HCPCS: E0148



- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 500 lbs
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl-contoured hand grip

PRODUCTS	Description
10220-1	Walker,Oversized,Adt.,1/ea

SPECIFICATIONS	Description
Carton Shipping Weight	12 21 18 lbs.
Overall Width	25.25" 24.5"
Warranty	Limited Lifetime
Weight	9 8 lbs.
Weight Capacity	500 lbs
Width	20.5" inside hand grip 23" inside back legs

ACCESSORIES	Description
10125HD	Walker Tray, HD, 1/EA, CS