

Nitro Elite CF, Carbon Fiber Rollator



HCPCS:




FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F




FIGURE G

- Luxurious super lightweight carbon fiber frame (Figure A)
- Frame weighs only 12 lbs for powerful performance with swift and effortless mobility
- Opulent and removable zippered storage bag Unique bag attachment keeps bag securely in place when rollator is open or folded (Figure B)
- Brake cable inside frame for added safety (Figure C)
- Handle height easily adjusts with unique push button
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort (Figure D)
- Easily fold with one hand to ultra-compact size for storage Handle is located on seat (Figure E)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable
- EZ clip cane holder included

PRODUCTS

Description

RTL10266CF	Nitro Carbon Fiber Rollator, Black
------------	------------------------------------

SPECIFICATIONS

Description

5 Year Limited Warranty	Brake and Handle Housing
Carton Shipping Weight	19 lbs.
Handle Height	33.25" - 38.2"
Note	*excluding bag.
Seat Dimensions	17.75" (W) x 20.75" (H) x 8.5" (D)
Unit Dimensions	23" (W) x 27.75" (L)
Warranty Limited Lifetime	Frame
Weight	12 lbs.*
Weight Capacity	300 lbs

PARTS

Description

1026638	Safety Lock, Nitro, 1/ea, CS
1026636	Step Tube, Left, Nitro, 1/ea, CS
1026637	Step Tube, Right, Nitro, 1/ea, CS
1026609-FSW-HEMI	SeatLiftingStrapW/AttachHrdwre,10266WT-, CS
1026643	Cross Brace Screw & Insert, Nitro, 1/ea, CS
1026644	Folding Lock Mechanism, Red, Nitro,1/ea, CS
1026641	WheelBrakeLeftW/Hrdwr,Plastic,Nitro,1ea, CS
1026618-AN	FrontWheelBearings,NitroW/NewFork, 1/pr, CS

PARTS	Description
1026642	WheelBrakeRightW/Hrdwr,Plastic,Nitro,1e, CS
10266CF-18	BrkHndlAsmby,10266CFBefrSer#2A1801000001, CS
10266CF-19	BrakeCableGuideAsmbly,Lft,RTL10266CF,1e, CS
10266CF-20	BrakeCableGuideAsmbly,Rgt,RTL10266CF,1e, CS
10266CF-22	Brake Cable, 10266CF, 1/ea, CS
10266CF-12	HandlebarAsmblyW/oBrake,Lt,RTL10266CF,1, CS
10266CF-13	HandlebarAsmblyW/oBrake,Rt,RTL10266CF,1, CS
10266CF-07	RearAxleBoltW/Nut&Washer,RTL10266CF,1/e, CS
10266CF-01	Back Assembly, RTL10266CF, 1/ea, CS
10266CF-04	FrontCasterW/Bearings,Left,RTL10266CF,1, CS
10266CF-03	Front Axle W/Hub Cap, RTL10266CF, 1/ea, CS
10266CF-05	FrontCasterW/Bearings,Rght,RTL10266CF,1, CS