

Long Term Care Bed

Model NV-LTCB-703



Features:

- Two way Red/Green color coded floor lock system
- Eight (8) Function Handheld
- Pendant Control
- Auto Contour
- 4 DC Motors Synchronized
- Hi/Lo control system keeps bed level
- Four sealed bearings secured in "G" track
- IV pole holders
- Extended Heel Lift
- Eight 3" Heavy Duty casters roll at any height
- Fixed adjustable wall bumper

Dimensions:

- Bed Height Range: 7" - 30"
- Bed Length Range: 76", 80" or 84"
- Bed Width: 36"
- Warranty: Frame 15 Year, Welds Lifetime, Electronics 3 Year
- Weight Capacity: 500 lbs
- Bed Height Min w/casters: 7"
- Bed Height Max w/casters: 30"

Options:

- Staff Control with Trend/Reverse Trend and Cardiac Chair
- 88" length Extension
- Deluxe Assist Rails, Assist Bars
- Six styles and twelve wood grain laminate colors for ends
- Battery Back-up, Bed Transporter, Trapeze