



NOVUM

Medical Products, Inc.

Long Term Care Adult Bed

Model NV-LTC-P503

Features:

- 450 lbs. weight capacity, Resident Weight 425 lbs
- Best in class height range 7.8" to 30"
- Integrated 76" and 80" length extension
- 8 Function Handheld Pendant Control includes Auto Contour
- 3 DC Motors keeps bed level with true vertical lift
- 4 sealed bearings for smooth/quiet lifting
- Eight 3" Heavy Duty casters roll at any height, 4 locking, 2 caster guide locks
- Fixed wall bumper

Specifications:

- Bed Width 36"
- Bed Height Range 7.8" - 30"
- Bed Length Range 76" or 80"
- Weight 210 lbs.
- Weight Capacity Accommodate up to 450 lbs. Maximum Resident Weight 425 lbs.
- Warranty 10 Years on Frame, 15 Years on Welds, 2 Years on Electronics

Options:

- Assist Rails, Assist Bars
- Bed ends, Six styles and twelve wood grain laminate colors



Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com