

## HDPE Maternity Bassinet

Model # 9201-C-A




### Features:


- Constructed of HDPE Plastic Resin for advanced infection control
- Easy to Clean
- Extremely Durable
- Easy to Maneuver
- Inhibits microbial growth
- Custom Dimensions available
- 5" Casters, two locking and two non-locking
- Quiet, simple reflux positioning
- Includes an 18.75" pull-out deck for extra work space and side chart holder

**Overall Dimensions:** 31"H x 33.5"W x 22"D


**Drawer Dimensions:** 9.5"W x 16"D




## HDPE Color Options



DOLPHIN GRAY




CHARCOAL GRAY



WHITE




BLACK



MOCHA BROWN




NUTMEG



INDIGO



SLATE GRAY



SANSHADE